

# Reflection Electron Energy Loss Spectroscopy UPGRADE for VersaProbe II, III and 4

**Reflection Electron Energy Loss Spectroscopy (REELS)** is compatible with *VersaProbe II*, *VersaProbe III*, and *VersaProbe 4*. REELS is a surface analysis technique in which a specimen is bombarded with an electron beam ( $\leq 2000$  eV) and the energy distribution of the reflected electrons is measured. This energy distribution contains features corresponding to discrete losses of energy of the reflected electrons due to excited atomic states, valence band transitions, material bandgaps, etc.


## Key Applications of REELS

- Analyze the electronics and bonding states of the materials
- Measure band gap of semiconductors
- Compare the relative concentration of hydrogen
- Discrimination of  $sp^2$  /  $sp^3$  character of carbon

## Analysis Examples


### Bonding Chemistry Analysis by REELS

REELS is sensitive to hybridization of carbon. Energy loss peak due to  $\pi \rightarrow \pi^*$  transition in graphite appears 6-8 eV lower than the elastic peak.




### Hydrogen Concentration Measurement by REELS

Hydrogen-derived energy loss peak is present due to the rebound effect of hydrogen atoms. Relative concentration of hydrogen, which is not accessible by XPS, can be obtained from the intensity of this energy loss.



## Band Gap Measurement by REELS

REELS spectrum of SiO<sub>2</sub> thermal oxide film (25 nm) on Si wafer. REELS can be used to calculate the band gap of SiO<sub>2</sub>. The band gap is obtained from the intersection of a straight line originating from the background level and a linear fit to the onset of loss signal spectrum as shown in the zoomed-in REELS spectrum.



## REELS Hardware

The REELS upgrade is a dedicated electron gun with a tungsten filament and is operated using *SmartSoft* software. REELS sensitivity on clean Au foil is  $\leq 0.5$  eV FAT (constant energy resolution) and  $\geq 1,500,000$  /nA can be operated up to 2 keV.

System requirements: Windows 10

