

# PHI 710 Scanning Auger Nanoprobe Consumable Parts List

---

## Contents

VACUUM SYSTEM .....	2
PUMPS.....	2
GAUGES.....	2
GASES.....	2
SEALS (Gaskets/O-Rings).....	2
COMPONENT PARTS .....	3
ANALYZER.....	3
ELECTRON GUN.....	3
SCINTILLATOR.....	3
ION GUNS.....	3
MISCELLANEOUS.....	4
FUSES.....	4
LAMPS/FOIL/BULBS.....	4
SAMPLE PLATES.....	4
COMPUTER/PRINTER .....	4

JULY 2017

# PHI 710 Scanning Auger Nanoprobe

Part Number	Description	Quantity
<b>VACUUM SYSTEM</b>		
<b>PUMPS</b>		
2141400	Sublimator Filaments	1
630660	Di Ion Pump Elements (2 upper pump, 6 lower pump)	8
727970(PHI) 320531(Gamma)	Triode Ion Pump Elements 240 l/s pump Note: MPC Negative polarity board (P/N 727917) required for changing from Di elements to Triode requiring 4ea retaining brackets (P/N 805422)	8
727971	Negative Polarity Board for MPC for triode pump elements	1
628239	Rough Pump Oil-Type P3-5 Liter	1 bottle
710201	Rough Pump Oil-(Ulvac Pump)-4 Liter	1 bottle
648914	Rough Pump-Oil Mist Filter-Copper Sieve, 4"	1
725038	Tip Seal Kit for Edwards Scroll Pump Model nXDS 15i	1
<b>GAUGES</b>		
705313	Convectron Gauge Tube, Varian	1
705319	Ion Gauge, Nude (Upper Nipple & Main Chamber)	2
725351	Filament Replacement-Thoriated-Iridium (Varian)	1
<b>GASES</b>		
627206	Argon Lecture Bottle-1400 PSIG	1
<b>SEALS (Gaskets/O-Rings)</b>		
725038	Tip Seal Kit for Edwards Scroll Pump Model nXDS 15i	1
715074	Seal Kit-VAT Gate Valve (V1 Seal Kit)	1
2680750	Gasket-Cu 1.33 OD Flange-Pkg 10	1
2681500	Gasket-Cu 2.75 OD Flange-Pkg 10	1
2682000	Gasket-Cu 3.38 OD Flange-Pkg 10	1
2682500	Gasket-Cu 4.50 OD Flange-Pkg 10	1

JULY 2017

# PHI 710 Scanning Auger Nanoprobe

2684000	Gasket-Cu 6.00 OD Flange-Pkg 10	1
2686000	Gasket-Cu 8.00 OD Flange-Pkg 10	1
D76219	Gasket-Cu 10.00 OD Flange-Pkg 1	1
627684	Gasket-Cu 4.50 OD/2.75 ID-Special gasket for Auger Analyzer	1

## COMPONENT PARTS

ANALYZER		
649372	Multi-Channel Detector-Channel Plates	1Set
628779	Multi-Channel Detector Anode-8 Channel	1
ELECTRON GUN		
649384	Thermal Field Emitter	1
648718	Aperture Plate	1
SCINTILLATOR		
710819	Secondary Electron Detector-Scintillator	1
ION GUNS		
647104	06-350 Long Life Ionizer	1

JULY 2017

# PHI 710 Scanning Auger Nanoprobe

MISCELLANEOUS		
	<b>FUSES</b>	
609040	Fuse, Pico, 1A, 125V, Radial	5
171300	Fuse-3AG, 0.50A, 250 V, Slo-Blo	5
171302	Fuse-3AG, 1.0A, 250 V, Slo-Blo	5
171303	Fuse-3AG, 1.50A, 250 V, Slo-Blo	5
601265	Fuse-3AG, 5.0A, 250 V, Slo-Blo	5
628240	Fuse-3AG, 20A, 250 V, Very Fast-Blo	5
628527	Fuse, AB, 10A, 250 V	5
628528	Fuse, AB, 15A, 250 V	5
1003791	Fuse, AB, 0.25A, 250 V	5
	<b>LAMPS/FOIL/BULBS</b>	
700256	Aluminum Foil-18" Wide x 250 Foot Roll	1
710744	Lamp-Halogen, JCR 12V-100W (Schott Light Source)	1
	<b>SAMPLE PLATES</b>	
713036	Model 191HC-Recessed Specimen Holder	1
716012	Model 193G-Faraday Cup-Augur	1
713034	Model 190HC-Flat Specimen Holder	1
716014	Model 194D-Beam Diameter-Augur	1
716010	Model 202A-30° Sample Holder	1
713038	Model 260HC-Sample Holder-Flat-0.6cm	1
710896	Holder-Sample, 710-HERO	1
	<b>COMPUTER/PRINTER</b>	
	ALL PC AND PRINTER MEDIA ARE CONSUMABLE	

JULY 2017