



18725 Lake Drive East,
Chanhassen, MN 55317
Phone: 952-828-6200
Fax: 952-828-6176
www.phi.com

January 13, 2022

Memo: Regarding Apache Log4j software library

Physical Electronics software packages, which include system control software, data reduction software, and standalone helper applications don't use the Apache's Log4j software library.

Please be assured that PHI's software will not be subjected to the vulnerability threat by the Apache Log4j library.

Please contact service@phi.com if you have any further questions or concerns.

Regards,

Melinda Caouette

Principal Software Engineer/Team Leader

Physical Electronics

18725 Lake Drive East, Chanhassen, MN 55317

