


## MultiPak 9.9.2

### New Features and Improvements


**MultiPak 9.9.2 for XPS and AES** is the latest update of our comprehensive data reduction and interpretation package.

The update includes new features and improvements:

- Φ New File Selection User Interface** contains tooltips information, minimized number clicks, more efficient control of folders hierarchy, adjustable column width, and more.




- Φ Optimized upper toolbar control** for most frequently used functionalities.




- Φ Updated display of header** – TOA and sputter time for spectral data, the field of view for SXI and MAPs.
- Φ Convenient control of windows** - minimization of all windows by minimizing the Spectral window.
- Φ Saving current cycle/TOA or line** as an individual SPE file from PRO/ANG/LIN files.
- Φ Improved export options.**
- Φ Kinetic energy display** on peak ID cursor.

Φ Open XPS and AES files in the same work session.



Φ Customizable markers on profile displays.

Φ Improved Smart background functionality



Φ Automated work function/ionization potential calculations from UPS spectra.

Φ PowerPoint™ connect is enabled by default.